

ARTHRAMID®


What is Arthramid®?

Arthramid® is a unique patented 2.5% injectable polyacrylamide hydrogel (2.5% iPAAG) that is injected directly into the joint to improve mobility and reduce discomfort.

Arthramid® goes beyond other therapies, incorporating a dynamic bioscaffold technology to safely and effectively manage osteoarthritis long-term.

How Does Arthramid® Work?

Arthramid® targets the joint's synovial membrane - the soft tissue that produces joint fluid. After injection, it gradually integrates into this membrane, a process that takes approximately 2 weeks.



Through this action, Arthramid® reinforces and supports the joint capsule to:

- Reduce stiffness and restore flexibility
- Enhance stability in the joint
- Improve joint fluid quality, which cushions, nourishes cartilage, and reduces friction
- Help break the cycle of arthritis and optimize overall joint health long-term

With Arthramid®, your dog can enjoy a more active, comfortable life - helping them get back to running, playing, and moving with ease.

Why Vets and Dog Owners Trust Arthramid®

- Long-lasting relief from osteoarthritis
- Improves mobility and joint health
- Safe, non-invasive, and backed by 20+ years of research
- Works for all stages of osteoarthritis, from early to advanced
- Reduces the need for daily pain medication or frequent injections
- Safe for use in multiple joints

PROVEN RESULTS: WHAT THE STUDIES SHOW ABOUT ARTHRAMID®

- **85% of dogs** show relief from the clinical signs of osteoarthritis¹
- **63% of dogs** reduce or discontinue the use of other pain medications¹
- **Improved mobility & range of motion** in treated joints
- **Long-lasting benefits** with increased freedom of movement and overall improvement in quality of life over time

Testimonials

“Arthramid® has transformed how I treat osteoarthritis. It’s easy to administer, highly effective, and has redefined my approach to post-op recovery and OA management.”

Dr. Julie Letosky, DVM, CCRP
Up and Running Canine Rehabilitation

“Arthramid® should be a go-to for earlier OA management in small animals. It’s safe, long-lasting, and supports healthier joints.”

Dr. Rachel Horsley, DVM
Dogwood Veterinary Hospital

Post-Treatment Care

Your dog's joint comfort will gradually improve as Arthramid® takes effect.

The gel integrates into the joint's synovial membrane within 2 weeks, but visible improvement in your dog's movement and comfort may take up to 6 weeks. Trust the process and allow time for the treatment to deliver its full benefits.

HOW TO SUPPORT YOUR DOG'S RECOVERY POST-INJECTION

To help your dog get the most out of their treatment, follow these simple steps:

Rest for 48 hours: Keep activity minimal for the first two days

Controlled exercise: Stick to short, gentle leash walks for the next 2 to 4 weeks, slowly increasing duration as your vet advises

Vet check-up: Schedule a re-examination with your vet 4 to 6 weeks after treatment

Continue rehab: If your dog has a rehab plan, keep following it as recommended

